## sharpsafe

**Never Overfill Your Sharps Container** 

Never fill above the fill line. Never place your hands inside the container. ✓ MEDICINAL SHARPS Never touch the base of the container, always use the handle. DANGER Do not shake the container or try to force contents down. Familarise yourself 5 with your local Sharps **Management Policies** and Risk Assessments.

## Safer today for a sustainable tomorrow.

1 Western Avenue, Matrix Park, Buckshaw Village, Chorley, PR7 7NB, UK +44 (0)1204 529494 hello@vernacare.com vernacare.com

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. © Copyright Vernacare International Limited 2024. All rights reserved. CWM036 REV00



## Vernacare

WARNING: DO NOT FILL Above the line

TO BE

clinical waste management